

Import Minimac Output

Author: Jesse Dupre, Golden Helix, Inc.

Overview

This script will import Minimac info and dose files with phased genotype dosages which are output from running the Minimac software. For more information on Minimac see:

<http://genome.sph.umich.edu/wiki/Minimac>

Recommended Directory Location

Save the script to the following directory:

***..\Application Data\Golden Helix SVS\UserScripts\SVS\Import**

Note: The **Application Data** folder is a hidden folder on Windows operating systems and its location varies between versions. The easiest way to locate this directory on your computer is to open SVS and select the **Tools >Open Folder > User Scripts Folder** menu option. If saved to the proper folder, this script will be accessible from the project navigator **Import** menu.

Using the Script

1. From an open project select **Import > Import Minimac Output**
2. Select the *.info file(s) and optionally the *.dosage file(s). If dosage files are imported, the markers can be filtered based on R² values.

Two spreadsheets are created, one for the info file and the other containing the dosage data.